

# Measuring Microscopes Hyper MF / MF-U

CATALOG No. E4267-176



Measuring Microscopes So Accurate They Defy  
Common Sense

**Mitutoyo**

# World's Highest Measuring Accuracy



## Concept

Inspecting complex microstructures of ever-decreasing size demands ever-higher accuracy from measuring microscopes used to satisfy the manufacturing and quality control principle of Observation plus Measurement. Mitutoyo is committed to providing microscopes that meet this requirement as well as exceeding users' expectations in terms of sophisticated functionality and ergonomic features that allow fatigue-free use over extended periods of time.




## Core Technology

Over many years Mitutoyo has made significant contributions to the technologies that are key to the core technology of manufacturing industry: measurement. The experience and expertise gained is reflected in the design and manufacture of each individual component of these microscopes and can be seen most clearly in their sublime integration of optics, mechanics, and electronics.




## Ergonomic Design

The microscope main unit has been designed with the emphasis on user friendliness and ease of operation. Mitutoyo has executed the mechanical design to allow easy operation. Even after extended use, its fatigue-fighting design still provides a comfortable work-experience for the operator.



## World's Highest Measuring Accuracy\*

Measuring accuracy in the X- and Y-axes at full stroke surpasses class 0 of the JIS Standard for measuring microscopes (B7153-1995). This makes these microscopes ideal for high accuracy measurement of precision molds or cutting tools that require the best resolution, or for inspecting sub-miniature semiconductor / electronic parts such as wafers and integrated circuits.



\*As of April, 2004

## Large, Highly Accurate XY Stage

Mitutoyo uses a type of linear guideway on the large XY stage that is highly regarded for excellent straightness and stability. This is one key element in the strategy to maximize geometrical accuracy - another is FEM analysis. Our designers used FEM techniques extensively during the design phase to ensure stage stability was optimal in any measurement situation. Thus, the foundations for achieving the highest measuring accuracy were laid.



FEM analysis of the XY stage

## More about the XY Stage

The XY stage is a massive, highly stable design created using mechanical techniques developed over Mitutoyo's long years of experience in manufacturing precision measuring microscopes. Maximum stage loading is 30kgf (66lbf) and a range of useful fixtures is available that includes a wafer holder and swivel-center support.



## Highly Accurate Digital Scales

These microscopes are equipped with highly accurate digital glass scales on all three axes. Mitutoyo produces glass scales in an underground laboratory where the temperature and humidity are constant throughout the year. The XY (stage) and Z (optical tube) displacements are displayed digitally.




# Excellent Operability and Solid Reliability



## FS Optical System

The FS optical system is respected more than ever before for its ability to enable measurement, observation and analysis with a leading-edge combination of long working distance and high NA. This optical system ensures high operability when measuring deep holes, steps, etc., or when setting up

workpieces for measurement.



**Mitutoyo**

## Tilting Optical Tube\*

To reduce fatigue due to extended use, it is important that the operator use a microscope in an unforced posture. The eyepiece unit allows stepless adjustment of tilt angle so that, no matter what their physique, operators can always adjust the viewing position for comfortable working during any measuring task.

\* Available for model MF-U only.



## LAF Optical Tube\*


The LAF (Laser AF) optical tube can be selected as an option. The LAF system achieves high repeatability when measuring minute steps, etc., enabling difficult measurements with minimum fatigue.

\* Available for model MF-U only



## Front Operation

Controls are arranged to fall within easy reach of the operator on two control panels at the front of the microscope. This allows the operator to concentrate on measurement without having to look away from the eyepieces. Membrane technology makes the switches very durable.




# Measuring with Hyper MF - an Emotional Experience



## Cross-hair Reticle

How accurately the reticle can be aligned with a workpiece feature is a very important feature in a measuring microscope. Taking ocular resolution into account, the thin-line reticle has been standardized on a broken, 90° cross hair with a line width of 5 $\mu$ m\*. This allows precise positioning of the reticle.



\* Cross-hair reticles of 3 $\mu$ m and 7 $\mu$ m line width are also available.

## Fiber-optic Cold Light Illumination

A fiber-optic cold light illuminator and an IR absorption filter greatly reduce thermal effects on the instrument and workpiece that would otherwise have an adverse effect on measuring accuracy. Telecentric illumination is used for reflected light observation and Koehler illumination for viewing contours. Both systems use an aperture diaphragm for even, glare-free illumination with good image contrast.



## System Extensibility

A video port is standard on the optical tube, thereby allowing a digital camera unit and various vision analysis units to be added to a system.



Digital Color Vision System DV-520  
(Microscope main unit: MF-UB1010TH)

# Main Specifications

MF/MF-U

## Highest-in-Class Measuring Accuracy

Mitutoyo has achieved a measuring accuracy of  $\pm (0.9+3L/1000)\mu\text{m}$  (L: Measured length in mm) in both X- and Y-axes. This performance surpasses Class 0\* of JIS B7153:1995, Measuring Microscopes, and allows support of ultra-precise inspection and measurement of the smallest visible features to those extending across the full measuring range of these microscopes.

\* Class 0:  $(2+0.01L)\mu\text{m}$  or less, L: Measured length (mm)



MF/MF-U

## Three-axis Motor-driven Joystick Ensures High Operability from High Speed to Ultra-Low speed


The X, Y, and Z axes are driven and controlled with one joystick that serves as the nerve center of front operation. Speed control is possible from high-speed traverse of the stage to ultra low-speed, minute positioning of a workpiece. Also, the lock mechanism is provided for each X, Y, and Z axis to support high-accuracy pitch measurement by single-axis displacement. The primary target is assumed to concentrate the operator on a workpiece.



MF/MF-U

## Large, Highly Accurate XY Stage Handles Wide-field, Heavy-weight Workpieces


The pressures for diversification and up sizing of workpieces are increasing in various industrial fields, such as semiconductors, electronics, precision automotive parts and tools. These microscopes not only have the accuracy for the smallest workpiece but also have the power to handle larger components such as lead frames, precision cutting tools and molds.



MF-U

## Remote-controlled Objective Magnification Change

The power turret in the optical tube is controlled with membrane switches on the left front panel. LEDs indicating each lens position on the upper part of the optical tube are linked to rotation of the turret so that the operator can see the current magnification at a glance.



Power Turret Drive Switches

Left control panel

# Laser Auto Focus

## LAF Optional Tube


The laser auto focus function provides high accuracy and high repeatability and brings significant advantage to the inspection of minute steps, multi-layer board detail, etc. A powerful function that helps avoid operator error and ensures high productivity.



LAF is available both in BF and BD optical tubes.

## Selectable LAF Functions

Providing a choice of the Just Focus (JF) mode that functions quickly at the current point of interest and the Tracking Focus (TF) mode that tracks the focusing position to retain sharp focus as the stage moves has improved measurement efficiency.



Visible  
Semiconductor  
Laser 690nm

Image graphic of AF optical path



ultra-precision  
small-sized gear 10X



Wafer 20X



Metal slope 20X



Metal marking 20X

## LAF Effective in the Smallest Area

An LAF spot diameter of  $\phi 1\mu\text{m}$  or less is achieved using an objective with a magnification of 50X or more. This performance supports a wide range of measurement tasks.


The spot diameters are a logical value determined by calculation.

Objective	Spot diameter
MplanApo 2X	$6\mu\text{m} / 630\mu\text{inch}$
MplanApo 5X	$6\mu\text{m} / 240\mu\text{inch}$
MplanApo 10X	$3\mu\text{m} / 120\mu\text{inch}$
MplanApo 20X	$1.5\mu\text{m} / 60\mu\text{inch}$
MplanApo 50X	$0.8\mu\text{m} / 3\mu\text{inch}$
MplanApo 100X	$0.6\mu\text{m} / 2\mu\text{inch}$

The AF function delivers highly repeatable focusing on areas with different surface textures and slopes.

## Laser Beam Class

The LAF (factory-fit option) function uses a low-power laser that corresponds to Class 2 (visible light) of JIS C6802/1997, Safety of Laser Products.



# Main Unit Specifications

Hyper MF



Hyper MF-U




Model name		Hyper MF	Hyper MF-U
Order No.		<b>176-401</b>	<b>176-402</b>
Optical tube		-	(Bright field / bright field LAF / bright and dark field / bright and dark field LAF)* <sup>1</sup>
Standard reticle (built-in)		90° broken cross-hair (line width: 5μm) Note: Not selectable for Hyper MF-U	
Pupil distance adjustment		Siedentoph type Adjustment range: 51 to 76mm / 2.01" to 2.99"	
Optical path transfer ratio		Observation / video camera = 50/50	
Angle from horizontal		25°	Manual tilting
Video port		Standard equipment	
Observation image		Erect image	
Eyepiece Magnification		10X, 15X, 20X	
		Monocular unit (with one eyepiece) or Binocular unit (with two eyepieces) selectable	Provided with two 10X eyepieces
Objective	ML objective	1X, 3X (standard accessories), 5X, 10X, 20X, 50X, 100X	-
	Bright field (BF)	-	MplanApo, MplanApo SL, G Plan Apo
	Bright / dark field (BD)	-	BD Plan Apo, BD Plan Apo SL
Turret	Bright field (BF)	-	(Centering, parfocal, manual control, 4-way / manual control, 4-way with sensor / centering, motor-driven control, 5-way with sensor* <sup>2</sup> )
	Bright / dark field (BD)	-	(Manual control, 4-way / manual control, 4-way with sensor / motor-driven, 4-way with sensor* <sup>3</sup> )
Focusing Unit	Maximum workpiece height	150mm / 5.91"	
	Measuring accuracy* <sup>4</sup>	±(1.5+10L/1000)μm L: Measured length in mm	
	Drive method	Motor-driven control with the joystick	
Illumination	Contour illumination	Koehler illumination, aperture diaphragm built-in, 12V50W halogen bulb, 100-step light adjustment, fiber-optic cold light illumination	
	Surface illumination	Telecentric illumination, variable aperture diaphragm mechanism, 12V100W halogen bulb, 100-step light adjustment, fiber-optic cold light illumination	
Stage	Measuring range (X-axis, Y-axis)	250mmx150mm / 9.84"x5.91"	
	Measuring accuracy* (at no load on X- and Y-axis)	±(0.9+3L/1000)μm L: Measured length in mm	
	Top face dimensions	460mmx350mm / 18.11"x13.78"	
	Effective dimensions of stage glass	300mmx200mm / 11.81"x7.87"	
	Swivel range	±3°	
	Maximum loading	30kgf / 66lbf	
Detector unit		Motor-driven control with the joystick	
Detector unit		High-accuracy digital scale (Patent registered)	
Digital display	Resolution	0.01μm / .0004inch	
	Display axes	X, Y and Z axis	
	Data processing unit	QM-DATA200 or Vision Unit	
	Joystick lock	✓	
Control unit	Fine pitch	✓	
	Data output	✓	
	Digital display reset	✓	
	Illumination brightness controller	✓	
	LAF (Just focus)	-	✓ (with LAF installed)
	LAF Tracking focus	-	✓ (with LAF installed)
Turret remote-control		-	✓ (with the power turret installed)
External dimensions	Microscope main unit	880mmx913mmx730mm / 34.65"x35.94"x28.74"	880mmx913mmx770mm / 34.65"x35.94"x30.31"
	Power unit	160mmx476mmx381mm / 6.30"x18.74"x15"	
Mass	Microscope main unit	Approx. 250kg / 551.2lb	Approx. 255kg / 562.2lb
	Power unit	14kg / 30.86lb	
Power supply		100, 110, 120, 220, 230, 240VAC 50/60Hz Maximum power consumption 700W	

\*1, 2, 3: Factory-fit option

\*4, 5: The measuring accuracy conforms to JIS B 7153.

Upon bulb replacement, order the contour illumination halogen bulb 12V50W (02APA527) or surface illumination halogen bulb 12V100W (517181). The high-brightness specification bulb (12BAD602) is also available as a 12V100W bulb.

# System Configuration




# Data Processing System - 2D Calculating System -



## Specifications

Order No.	<b>264-149</b>
Display languages	Japanese, English, German, Italian, Spanish, Portuguese, Czech, Chinese, and Korean
Measurement value unit	Length: mm, angle: degree / degree, minute, second (selectable)
Resolution	0.01μm
Programming function	Create / run / edit measuring procedures
Statistical processing	Measurement item, quantity of data, maximum value, minimum value, mean, standard deviation, range, histogram
Number of elements in memory	Maximum 1000 elements
Element call	Point, line, circle, distance, ellipse, square hole, slotted hole, point and angle of intersection
Element key-in	Point element line element, circle element
Display unit	Monochrome graphic LCD (320x240 dots, backlighted)
Input connector	(1) RS-232C [2]: Connected to the counter on the main unit. (2) X, Y, Z: Linear scale input (for the types with no counter) (3) OPTOEYE: OPTOEYE edge signal (for connecting OPTOEYE M2) Note: Available only if linear scale input is used. (4) FS: Connected to the foot switch.
Output connector	(1) RS-232C [1]: Connected to an external PC (measurement results) (2) PRINTER: Connected to the printer / external printer (measurement results) (3) FD: Connected to the FDD unit (measurement result file, measuring procedure file)
Measurement result file output	RS-232C output (CSV format, MUX-10 format)
Power supply	100 to 240VAC, 50/60Hz (AC adapter used)
Maximum power consumption	24W (excluding optional accessories)
External dimensions	Approx. 250x242x310mm (including the stand) / 10.24"x9.53"x12.20"
Mass	Approx. 2.2kg / 4.85lbs

## Dimensions



## Features

- > Powerful 2D measurement capabilities with graphic display functions that make the most of the large LCD screen
  - > Graphical help on the screen guides the operator during measurement sequences.
  - > Measurement results are displayed automatically
  - > Measurement procedures (Part Programs) can be learnt by the system and easily repeated with position navigation help on screen
  - > Frequently-used combination measurements (e.g. circle-to-circle) are single-key operations
  - > The Automatic Identification (AI) function recognizes the feature type automatically, making preselection unnecessary
  - > Macros to initiate learned measuring sequences can be created at a keystroke
  - > Custom menus to suit specific requirements can be created
  - > Tolerance comparisons and various statistical evaluation options are possible for every measurement result.
  - > Measurement results can be output to MS Excel®\* in table form (CSV)
  - > Measurement results and Part Programs can be stored on the USB FDD unit or USB-Memory stick available as an optional accessory
  - > A free-standing table version with tilting device is available
  - > The next measurement can be started even while the last is printing out
- \* MS-Excel® is a registered trademark of Microsoft Corporation.

## USB-Memory Stick

### No. 12AAH034

- Used for saving / reading files of part programs, user macros, measurement results, etc.



## USB-FDD Unit

### No. 12AAH035

- Used for saving / reading files of part programs, user macros, measurement results, etc.



## Printer

### No. 12AAD032

- Used to print measurement results.



### [Optional accessories]

#### No. 908353

Printer paper (5 rolls)


\* An external printer (color/monochrome) compatible with ESC/P is also available.

Printer control code system: ESC/P, compatible with MS-DOS.

Pin-out: 24 pins

ESC/P printer cable (No. 12AAA804): 2m - - Option

Receipt printer print example



Print method	Serial-matrix thermosensitive method
Number of print digits	40 digits
Print speed	Maximum 52.5cps (normal character)
External dimensions (WxDxH)	160x170x65.5mm (printer main unit) / 6.30"x5.59"x2.58"
Standard accessories	Printer cable, printer paper (1 roll), AC adapter (for 100V)




## Features


- > Automatic edge detection tool and measurement macro icons enable single-key measurement
- > Graphics and measurement navigation function support ease of use
- > Image capturing / saving function
- > Tolerancing calculation results and statistical processing for each item
- > Measurement results can be exported to MS-Excel®\* in CSV format (allows unique inspection sheets to be created on the same PC)
- > Supports total measurement on a single screen.
- > Automatic light equalizing function faithfully reproduces illumination conditions

\* MS-Excel® is a registered trademark of Microsoft Corporation.

## Measuring a Workpiece Feature



1) Display the feature to be measured on the monitor, adjust the illumination and focus with the microscope main unit, and then select the feature and the edge detection tool.



2) Click in the vicinity of the feature with the mouse to automatically detect its edge and perform the measurement / calculation.



3) The measurement results for the feature selected are displayed on the monitor.

## Specifications

Order No.	<b>359-680</b>
Optical system	0.5X (0.5X video adapter built-in) when Mitutoyo measuring microscope is installed
Image detecting unit	High-sensitive 1/3" CCD color camera
External dimensions of CCD camera	160x170x65.5mm (printer main unit) / 6.30"x5.59"x2.58"
Mass of CCD camera	0.4kg / .88lb
Monitor magnification / monitor visual field (at the use of 3X lens)	Approx. 63X / 3.1x2.3mm / .48"X / .12"x.09"
Resolution	0.01μm / .4μinch
Measuring accuracy	Depends on the accuracy of the measuring microscope.
PC	
CPU	Pentium 4 2.26GHz 256MB memory
Display	15 inch LCD
OS	MS-Windows XP*
Software	QSPAK
Maximum power consumption	Max. 400W

\* MS-Windows XP and MS-Excel® are registered trademarks of Microsoft Corporation. The PC configuration and specifications are subject to change without notice.



# Lens and Illumination

## Eyepieces

### Monocular unit MF



**No. 176-302**

With one eyepiece 10X / 24

### Binocular Unit MF



**No. 176-303-1**

With two eyepiece 10X / 24

### Eyepiece MF/MF-U



Part name	WF10X / 24	WF15X / 16	WF20X / 12
No. (1-piece pack)	<b>378-856-5</b>	<b>378-857-5</b>	<b>378-858-5</b>
No. (2-piece pack)	<b>378-856</b>	<b>378-857</b>	<b>378-858</b>
Magnification	10X	15X	20X
Field number	24	16	12


### Double-image eyepiece MF



**No. 375-044**

Aids accurate measurement of hole-to-hole distances, hole diameters and sections using the double image generate when the feature under inspection is not aligned with the optical axis of the microscope.

- Magnification: 10X
- Field number: 22



### Rotary template eyepiece MF



**No. 176-357**

Uses a template to superimpose screw thread forms (metric UST, Uni-field Screw Thread), concentric circles, and crosshairs on the workpiece image for quick and efficient measurement by comparison.

- Magnification: 10X
- Field number: 21

### Protractor eyepiece MF



**No. 375-043**

Measures angle between workpiece edges by successive alignment with two crosshairs whose separation is adjustable and calibrated through 360°

- Field number: 21
- Resolution: 5'



ISO metric screw pitch (for 3X objective)	0.25, 0.3, 0.4, 0.45, 0.5, 0.6, 0.7, 0.75, 0.8, 0.9, 1, 1.25, 1.5, 1.75, 2, 2.5, 3, 3.5, 4
ISO unified screw thread / inch (for 3X objective)	80, 75, 64, 56, 48, 44, 40, 36, 32, 28, 24, 20, 18, 16, 147, 13, 12, 11, 10, 9, 8, 7, 6
Concentric circles with cross hairs (for 3X objective)	ø1, ø2, ø3, ø4, ø5
Dotted-line cross scale	60°, 90°

### Digital protractor eyepiece MF

**No. 176-313**

Measures angle between workpiece edges by successive alignment with reticle cross hairs whose rotation is digitally calibrated. Switching or resetting the resolutions is controlled with the standard accessory counter. Data output to an RS-232C equipped PC is possible.

- Magnification: 10X
- Field number: 18
- Reticle: 90° solid line, 45° broken line
- Angular resolution: 0.00° or 1°
- Power supply: 9VAC, 600mA
- Maximum power consumption: 4W
- Maximum angle value: ±369.99° or ±369.59°



## Objective

### ML Objective

MF

Model name	ML1X	ML3X	ML5X	ML10X	ML20X	ML50X	ML100X
Order No.	<b>375-036-1</b>	<b>375-037-1</b>	<b>375-034-1</b>	<b>375-039</b>	<b>375-051</b>	<b>375-052</b>	<b>375-053</b>
Magnification	1X	3X	5X	10X	20X	50X	100X
Numerical aperture N.A.	0.03	0.07	0.13	0.185	0.42	0.55	0.7
Working distance WD	61mm / 2.40"	77mm / 3.03"	61mm / 2.40"	51mm / 2.32"	20mm / .79"	13mm / .51"	6mm / .24"
Focal depth	306µm / 12.05µinch	34µm / 1.33µinch	16.3µm / .64µinch	6.2µm / .24µinch	1.6µm / .06µinch	0.9µm / .94µinch	0.6µm / .02µinch

### FS Objective Bright field (BF)

MF/MF-U

Model name	MplanApo 1X	MplanApo 2X	MplanApo 5X	MplanApo10X	MplanApo 20X	MplanApo 50X	MplanApoHR50X	MplanApo 100X	MplanApoHR100X
Order No.	<b>378-800-3</b>	<b>378-801-3</b>	<b>378-802-6</b>	<b>378-803-3</b>	<b>378-804-3</b>	<b>378-805-3</b>	<b>378-814-4</b>	<b>378-806-3</b>	<b>375-815-4</b>
Magnification	1X	2X	5X	10X	20X	50X	50X	100X	100X
Numerical aperture N.A.	0.025	0.055	0.14	0.28	0.42	0.55	0.75	0.7	0.9
Working distance WD	11mm / .43"	34mm / 1.33"	34mm / 1.33"	34mm / 1.34"	20mm / .79"	13mm / .51"	5.2mm / .20"	6mm / .24"	1.3mm / .05"
Focal depth	550µm / 17.32µinch	91µm / 3.58µinch	14µm / .55µinch	3.5µm / .14µinch	1.6µm / .06µinch	0.9µm / .04µinch	0.48µm / .02µinch	0.6µm / .02µinch	0.34µm / .01µinch
Model name	MplanApoSL20X	MplanApoSL50X	MplanApoSL80X	MplanApoSL100X	MplanApoSL200X	GplanApo20X(t2.5)	GplanApo50X(t3.5)		
Order No.	<b>378-810-3</b>	<b>378-811-3</b>	<b>378-812-3</b>	<b>378-813-3</b>	<b>378-816-3</b>	<b>378-847</b>	<b>378-848-5</b>		
Magnification	20X	50X	80X	100C	200X	20X	50X		
Numerical aperture N.A.	0.28	0.42	0.5	0.55	0.62	0.28	0.5		
Working distance WD	30.5mm / 1.2"	20.5mm / .81"	15mm / .59"	13mm / .51"	13mm / .51"	29.42mm / 1.16"	13.89mm / .55"		
Focal depth	3.5µm / .14µinch	1.6µm / .09µinch	1.1µm / .04µinch	0.9µm / .04µinch	0.7µm / .03µinch	3.5µm / .14µinch	1.1µm / .04µinch		

### FS Objective Bright / dark field (BD)

MF/MF-U

Model name	BDplanApo2X	BDplanApo5X	BDplanApo10X	BDplanApo20X	BDplanApo50X	BDplanApoHR50X	BDplanApo100X	BDplanApoHR100X
Order No.	<b>378-831-5</b>	<b>378-832-6</b>	<b>378-833-5</b>	<b>378-834-5</b>	<b>378-835-5</b>	<b>378-845</b>	<b>378-836-5</b>	<b>378-846</b>
Magnification	2X	5X	10X	20X	50X	50X	100X	100X
Numerical aperture N.A.	0.055	0.14	0.28	0.42	0.55	0.75	0.7	0.9
Working distance WD	34mm / 1.34"	34mm / 1.34"	34mm / 1.33"	20mm / .79"	13mm / .51"	5.2mm / .20"	6mm / .24"	1.3mm / .05"
Focal depth	91µm / 3.58µinch	14µm / .55µinch	3.5µm / .15µinch	1.6µm / .96µinch	0.9µm / .04µinch	0.48µm / .03µinch	0.6µm / .02µinch	0.34µm / .01µinch
Model name	BDplanApoSL20X	BDplanApoSL50X	BDplanApoSL80X	BDplanApoSL100X				
Order No.	<b>378-840-5</b>	<b>378-841-5</b>	<b>378-842-5</b>	<b>378-843-5</b>				
Magnification	20X	50X	80X	100C				
Numerical aperture N.A.	0.28	0.42	0.5	0.55				
Working distance WD	30.5mm / 1.2"	20mm / .79"	13mm / .51"	13mm / .51"				
Focal depth	3.5µm / .14µinch	1.6µm / .09µinch	1.1µm / .04µinch	0.9µm / .04µinch				

### Twin fiber-optics illuminator

MF



No. 176-416

This uses the surface illumination light source in microscope main unit. The light equalizing function and condenser lens are included. 12V100W

\*10X lens or less is applicable.

### Fiber-optics ring-light

MF/MF-U



No. 172-417

This uses the surface illumination light source in microscope main unit. The light equalizing function and condenser lens are included. 12V100W

\*10X lens or less is applicable.

### LED ring-light

MF



No. 172-367  
(white LED)

No. 176-368  
(Green LED)

Position is adjustable so as to be appropriate for the light equalizing function and working distance. 12V 7.7W, outside diameter: 70mm/2.76"

\*10X lens or less is applicable.

### LED ring-light

(for FS Objectives)

MF-U



No. 78829  
(white LED)

Position is adjustable so as to be appropriate for the light equalizing function and working distance. 12V7.7W, outside diameter: 70mm/2.76"

\*10X lens or less is applicable.

# Main Optional Accessories

## Digital Color Vision



C-mount compatible high-definition digital color vision system. Allows anyone to easily capture and save high quality images from the microscope.

Three types: 5.2M pixel, 2.1M pixel, and 0.5M pixel.

•Main functions Image filing, on-screen manual measurement, white balance, comments on an image, scale display, etc.


**No.176-380:** DV-50CL (LCD spec)

**No.176-384:** DV-210L (LCD spec)

**No.176-382:** DV-520L (LCD spec)

\* On-screen manual measurement requires the calibration chart (**02AKN020**) separately.

## 3D Vision Combination Program VCPAK



This system is based on an omnifocal image focused in the whole visual field of the microscope to greatly improve observation at small depths of focus. It also enables 3D vision display.

•Main functions High-quality digital imaging, real-time vision combination by capturing a maximum of 256 sheets of images, adoption of universal image file formats (JPEG, BMP, TIFF, DIB), thumbnail display of saved images, noise reduction function, etc.

## Polarization Unit



**No.378-092**

Bright field (BF)  
Bright and dark field (BD)

## DIC Unit



**No.378-080:** For 5X, 10X

**No.378-079:** For 20X

**No.378-078:** For 50X, SL20X

**No.378-076:** For 100X, SL80X, SL50X

## Illumination Filter



**No.12AAA643:**ND2

**No.12AAA644:**ND8

**No.12AAA645:**GIF

**No.12AAA646:**LB80

## Centering, Parfocal Turret



**No.378-018**

Objectives: Up to 4

Visual field adjustment range:

±0.5mm/.02"

Parfocal adjustment range:

±0.5mm/.02"

## Power Turret (for BF)



**No.176-411**

Objectives: Up to 5

Visual field adjustment range:

±0.5mm/.02"

Positioning accuracy:  $2\sigma=3\mu\text{m}/120\mu\text{inch}$

Drive life: 1,000,000 movements

## V Block with Clamp



**No.172-378**

Maximum clamp diameter:

ø25mm/.98"

Mass: 0.8kg/ 1.76lb

\* Stage Adapter A is used together.

## Swivel Center Support



A tilt angle of  $\pm 10^\circ$  can be supported.

Minimum reading of angle:  $1^\circ$

Maximum support size: ø80x140mm/

3.15"x5.51" in horizontal orientation

ø65x140mm/ 2.56"x5.52" at a tilt

angle of  $10^\circ$

Mass: 2.5kg/ 5.51lb

\* Stage Adapter A is required.

## Holder with Clamp



**No.176-107**


Maximum workpiece thickness:

35mm/ 1.38"

Mass: 0.4kg/ .88lb

\* Stage Adapter A is required.

## Wafer Holder



**No.176-414**

Wafer size: 3 to 8 inches

Rotary Holder (**No.378-363**) attached

Mass: 3.2kg/ 7.06lb

## Center Support



**No.176-415**

Maximum support length:

250mm/9.84"

Maximum support diameter:

ø150mm/5.91"

Maximum diameter allowing external

apex observation: ø140mm/5.51"

Effective stroke of the center-clamping

mechanism: 22mm/ .87"

Mass: 13kg/ 28.66lb

## Stage Adapter A



**No.176-3042**

pieces pack Mass: 1.5kg/ 3.31lb

## Rotary table with fine feed knob (B)



**No.3176-306**

Effective glass diameter: ø240mm/9.45"


Rotary table rotation angle: Approx.

$51.5^\circ$  (per full turn of the fine feed

knob)

Mass: 6.5kg/ 14.33lb

## Vibration Damping Pad



**No.176-419**

Spring pad type

Float unit material: SUS304

WxDxH = 800x900x49mm/

31.50"x35.43"x1.93"

Mass: 58kg/ 127.9lb

## Machine Stand



**No.176-418**

Max. loading: 400kgf/ 881.8lbf

WxDxH = 1100x900x650mm/

433.07"x35.43"x26.59"

Mass: 45kg/99.21lb

# Measuring Microscope Line-up

Demand for measuring microscopes that can perform observational tasks as well as measurement is increasing rapidly in various sectors of industry such as semiconductors, electronic parts, precision auto parts and tools. The following summarizes Mitutoyo's line-up of measuring microscopes actively participating in many industries. Mitutoyo intends to widen the appeal of measuring microscopes that can determine miniscule part dimensions on a workpiece and make them serve as the Basic Machine for non-contact measurement.

## High-accuracy Measuring Microscopes

### Hyper MF/MF-U



Series name	Hyper MF	Hyper MF-U
Optical tube	Standard (finite distance compensation)	Metallurgical microscope (finite distance compensation)
Measuring range (X·Y·Z)	250·150·150mm	
Control system/ reading unit	3-axis motor-driven with Joystick/digital scale	
Resolution	0.01μm	
Data processing unit	QM-DATA200/vision unit	
Video port	Standard equipment	

## Measuring Microscopes

### MF-B/UB



Series name	MF	MF-U
Optical tube	Standard (finite distance compensation)	Metallurgical microscope (finite distance compensation)
Measuring range (X·Y·Z)	50·50·150/100·100·150/200·100·150/200·170·220/ 300·170·220/400·200·220mm	
Control system/ reading unit	Manual/digital scale	
Resolution	0.1 / 0.5 / 1μm	
Data processing unit	QM-DATA200/vision unit	
Video port	Standard equipment	

## Toolmaker's Microscope

### TM



Series name	TM
Optical tube	Standard (finite distance compensation)
Measuring range (X·Y·Z)	50·50·115/100·50·107mm
Control system/ reading unit	Manual/micrometer head
Resolution	1 (MHD head)
Data processing unit	QM-DATA200
Video port	None



Our products are classified as regulated items under Japanese Foreign Exchange and Foreign Trade Law. Please consult us in advance if you wish to export our products to any other country. If the purchased product is exported, even though it is not a regulated item (Catch-All controls item), the customer service available for that product may be affected. If you have any questions, please consult your local Mitutoyo sales office.

- Coordinate Measuring Machines
- Vision Measuring Systems
- Form Measurement
- Optical Measuring
- Sensor Systems
- Test Equipment and Seismometers
- Digital Scale and DRO Systems
- Small Tool Instruments and Data Management

Note: All information regarding our products, and in particular the illustrations, drawings, dimensional and performance data contained in this pamphlet, as well as other technical data are to be regarded as approximate average values. We therefore reserve the right to make changes to the corresponding designs, dimensions and weights. The stated standards, similar technical regulations, descriptions and illustrations of the products were valid at the time of printing. Only quotations submitted by ourselves may be regarded as definitive.

**Our products are classified as regulated items under Japanese Foreign Exchange and Foreign Trade Law. Please consult us in advance if you wish to export our products to any other country. If the purchased product is exported, even though it is not a regulated item (Catch-All controls item), the customer service available for that product may be affected. If you have any questions, please consult your local Mitutoyo sales office.**



**Mitutoyo Corporation**  
20-1, Sakado 1-Chome,  
Takatsu-ku, Kawasaki-shi,  
Kanagawa 213-8533, Japan  
T +81 (0) 44 813-8230  
F +81 (0) 44 813-8231  
<http://www.mitutoyo.co.jp>

